

Insio Nx

Technical Data

Made for É iPhone	iPad iPod		7Nx	5Nx	3Nx	2Nx	1Nx
ITE	ITC	CIC			IIC		

Insio Nx | Features and Accessories

	ITE / ITC				
	7Nx	5Nx	3Nx	2Nx	1Nx
Audiology					
Own Voice Processing (OVP) 1)	_	_	_	_	_
3D Classifier				_	_
Signal processing (channels) / Gain/MPO (handles)	48 / 20	32 / 16	24 / 12	16 / 8	16 / 8
Hearing programs	6	6	6	4	4
Sound Clarity					
HD Spatial	•	•	•	_	_
Extended dynamic range	•	•	•	•	•
Extended bandwidth	•	_	_	_	_
EchoShield	•	_	_	_	_
HD Music (presets)	3	1	_	_	_
eWindScreen binaural 1) 2)	•	•	_	_	_
eWindScreen	•	•	•	•	_
Noise Management					
Speech and noise management (steps)	7	5	3	3	1
SoundSmoothing (steps)	3	3	1	1	_
Directional speech enhancement (steps)	3	1	_	<u> </u>	_
Feedback cancellation	•	•	•	•	•
Speech Quality					
Directionality (channels) 6)	48	32	24	16	16
Binaural OneMic Directionality ^{1) 5)}	0	0	0		_
Narrow Directionality 1)	•	•	•	_	_
Spatial SpeechFocus 1) 3)	•	•			<u> </u>
SpeechFocus	•	•	_	_	_
TwinPhone ¹⁾	•	•	•	_	_
Frequency compression	•	•	•	•	•
Direct Streaming					
Made for iPhone	0		0	0	0
Adaptive Streaming Volume 4)		0	0		0
Tinnitus	O			0	
Notched Noise Therapy	•	•	•		_
Tinnitus noiser				•	_
Fitting					
Smart Optimizer and Data Logging	•	•	•	•	
Acclimatization manager	•	•	•	•	•
Performance Guide	•	•	•	•	•
Insitugram	•	•		•	•
Learning (classes)	6	3	1		_
TeleCare	-				
Basic Remote Tuning	•	•	•	•	•
Full Live Remote Tuning	\bigcirc	0	0	0	0
			_		

Insio Nx | Features and Accessories

			CIC / IIC		
	7Nx	5Nx	3Nx	2Nx	1Nx
Audiology					
Own Voice Processing (OVP) 1)	_	_	_	_	_
3D Classifier		_	_	_	_
Signal processing (channels) / Gain/MPO (handles)	48 / 20	32 / 16	24 / 12	16 / 8	16 / 8
Hearing programs	6	6	6	4	4
Sound Clarity					
HD Spatial	_	_	_	_	_
Extended dynamic range	•	•	•	•	•
Extended bandwidth	•	_	_	_	_
EchoShield	•	_	_	_	_
HD Music (presets)	3	1	_	_	_
eWindScreen binaural 1) 2)	_	_	_	_	_
eWindScreen	•	•	•	•	_
Noise Management					
Speech and noise management (steps)	7	5	3	3	1
SoundSmoothing (steps)	3	3	1	1	<u>—</u>
Directional speech enhancement (steps)	1	1	_	_	_
Feedback cancellation	•	•	•	•	•
Speech Quality					
Directionality (channels) 6)					
Binaural OneMic Directionality ^{1) 5)}	•	•	•	_	_
Narrow Directionality 1)	_	<u>—</u>	_	<u>—</u>	<u>—</u>
Spatial SpeechFocus 1) 3)	_	<u> </u>	_	—	—
SpeechFocus					
TwinPhone ¹⁾	•				
Frequency compression					
Direct Streaming					
Made for iPhone	_	_			
Adaptive Streaming Volume 4)		_			
Tinnitus					
Notched Noise Therapy	•	•		_	_
Tinnitus noiser	•	•		•	
Fitting					
Smart Optimizer and Data Logging	•	•	•	•	•
Acclimatization manager	•	•	•	•	•
Performance Guide	•	•	•	•	•
Insitugram	•	•	•	•	•
Learning (classes)	6	3	1	<u> </u>	_
TeleCare					
Basic Remote Tuning	•	•	•	•	•
Full Live Remote Tuning		_			

Insio Nx | Features and Accessories

	ITE / ITC (with BT)	ITE / ITC (without BT)	CIC	IIC		
	7Nx / 5Nx / 3Nx / 2Nx / 1Nx					
Style Specific Features						
Ingress Protection Rating	_	_	_	_		
Telecoil	_	0	_	_		
Charging contacts	_	_	_	_		
Battery Size	312	312	10	10		
Battery door on/off function	•	•	•	•		
Nanocoated housing		_	_	_		
e2e wireless 3.0	•	0	0	•		
User controls coupling via e2e	•	•	•	•		
Wireless programming	•	0	0	•		
Instrument configurations						
Flat cover		_	_	_		
Rotary volume control	_	0	_	_		
Push button	•	0	0	_		
Rocker switch	_	_	_	_		
Color conversion kit	_	_	_	_		
Battery door – direct audio input	_	_	_	_		
Battery door – child lock	_	_	_	_		
Small earhook	_	_	_	_		

Accessories				
miniPocket	\circ	0	0	0
CROS Silk Nx 7)	_	0	0	0
CROS Pure 312 Nx 7)	0	_	_	_
StreamLine TV	0	_	_	_
StreamLine Mic	0	_	_	_
Арр				
Signia App	\circ	0	0	0

¹⁾ req. bilateral fitting

■ available ○ optional — not available■ ■ ■ ■ highest feature performance

²⁾ not available in the universal program on 5Nx

³⁾ for 5Nx in Stroll Program or with Spatial Configurator only

⁴⁾ streaming only

⁵⁾ req. single mic instruments

⁶⁾ req. directional microphone

 $^{^{7)}}$ for 7Nx / 5Nx / 3Nx available only